


COVID-19 Disinfection for Eye Protection 11.3.20 1300

Spectrum Health contact: Jen Grile *Highlight denotes new content

Document reviewed: 3.25.2021

Seq. No	Task Description:	Key Point / Image / Measure (what good looks like)	Who
<p><i>The following steps should be followed for disinfection goggles, face shields, surgical eyewear, and safety glasses ONLY. Masks with face shields do not apply – these should be disinfected and stored following the mask guidelines.</i></p>			
1.	Obtain eye protection	<p>Goggles, face shield, or surgical eye wear</p>  <p>Face shields or goggles are required for the care of confirmed or suspected COVID-19 patients</p>	Health Care Worker
2.	Inspect eye protection for signs of deterioration (cracking, clouding, loss of sheen, compromised integrity, etc)	<ul style="list-style-type: none"> If found, discard eye protection 	Health Care Worker
3.	Perform hand hygiene prior to donning eye protection	<ul style="list-style-type: none"> Do not touch the front of the goggles, face shield, surgical eye wear, or safety glasses 	Health Care Worker
4.	Wear eye protection continuously and for multiple patient encounters	<ul style="list-style-type: none"> Follow extended-use guidelines 	Health Care Worker
5.	Disinfect eye protection immediately upon doffing	<ul style="list-style-type: none"> Don gloves Wipe the eye protection with hospital-approved disinfectant. <ul style="list-style-type: none"> Allow for contact time per label. 	Health Care Worker

Seq. No	Task Description:	Key Point / Image / Measure (what good looks like)	Who
		<ul style="list-style-type: none">○ Wipe off any residue with a paper towel.• Remove gloves.• Perform hand hygiene.• Place goggles, face shield, surgical eyewear, or safety glasses in isolation cart or on wall hook.	